


2. Steps of DNA cloning





2.1 Preparation of the plasmid

Plasmid minipreparation (质粒的小量制备) or miniprep:

Isolation of plasmid DNA from a few mililiters (mL) of bacterial culture.


CsCl gradient purification

Grow the cell Harvest the cell by centrifugation Resuspend (重悬) the cell pellet Alkaline lysis (碱裂解) of the cells Neutralization (中和) Phenol extraction (酚抽提) Ethanol precipitation(乙醇沉淀)




2.1.1 Resuspension

- Solution I 50 mM glucose, 25 mM Tris-HCl (pH 8.0), 10 mM EDTA (pH 8.0)
- Glucose prevents *E. coli* from quickly falling to the bottom of the tube, and the Tris-HCl buffers the cell at pH=8.0.
- After cell lysis (next steps), EDTA limits
 DNA degradation by binding (chelating, 整
 合) Mg²⁺ ions that are a needed cofactor for bacterial nucleases




2.1.2 Alkaline lysis

- Solution II 0.2 M NaOH, 1% SDS
- NaOH denatures DNA and hydrolyzes RNA.
- SDS disrupts cell membrane and denatures proteins. (碱性条件下裂解 效率高)




2.1.3 Neutralization

- Solution III KAc-HAc buffer (pH 5),
 5 M KAc
- KAc renatures plasmid DNA (supercoiled) and precipitates denatured proteins, along with the chromosomal DNA and most of the detergent (potassium dodecyl sulfate is insoluble in water).




2.1.4 Phenol extraction

Phenol extraction to get rid of the protein contaminants






2.1.5 Cesium chloride gradient purification




Supercoiled DNA may be purified from protein, RNA, chromosomal DNA and nicked plasmid DNA in one step!




2.1.6 Ethanol precipitation

- 0.3 M Na⁺(NaAc, pH 5.2), 3 volumes of ethanol
 - 乙醇消除核酸的水化层,使带负电荷的磷酸基团暴露出来。
 - Na+与磷酸基团结合,在沉淀形成部位降低多核苷酸链之间的排斥力。



- Centrifuge 13 000g and air dried
- Tris-EDTA (TE) buffer (EDTA chelates Mg²⁺)
- RNase A digests remaining RNA contamination



2.2 Get the target gene

(RT-)PCR amplification of the target gene

Chemical synthesis of the target gene


Restriction digestion from DNA library

Genomic library (基因组文库)



2.2.1 PCR amplification of the target gene










(2) Features of PCR

- 引物设计: 引物的长度一般为18~25 bp,且 引物自身及两个引物之间不能过多配对。3° 端不得有任何修饰。5°端加酶切位点。
- ▶ 扩增最有效长度:一般在200~800bp之间
- ▶ 扩増效率: 2ⁿ, 20次后不成此关系,因此 n=25~30次(荧光定量PCR n=40)。
- <u>变性温度与时间</u>: 94℃,第一次模板链较长,在5~10 min。每次循环中变性1~2 min。
- 退火温度:50~60℃,引物G+C含量一般在 40~60%之间,两条Tm值相差2~3℃以内。












(4) Detection of PCR products

琼脂糖凝胶电泳 (agarose gel electrophoresis)


- Agarose
 - Agarose is a polysaccharide derived from seaweed, which forms a solid gel when dissolved in aqueous solution at concentrations between 0.5% and 2% (w/v).
 - Agarose is a more purified form of the agar (琼脂) used to make bacterial culture plates.



Principle

Migration

When an electric field is applied to an agarose gel in the presence of a buffer solution, DNA fragments move through the gel towards the + electrode, because DNA negatively charged.



> Migration rate

The Migration rate is dependent on fragment sizes and shapes of the DNA and the electrophoresis may be used to separate mixtures of DNA fragments by the rate.

超螺旋>线型>开环

琼脂糖凝胶浓度	线性DNA的有效分离范围
0.20/	5 60 1-1-

0.5% 5-60 Kb 0.6% 1-20 kb

70/ 0.0.10.11


0.7% 0.8-10 kb

0.9% 0.5-7 kb


1.2% 0.4-6 kb


1.5% 0.2-4 kb

2.0% 0.1-3 kb




Staining for gel






The DNA is stained by the ethidium bromide (EB) in the gel, the DNA shows up as an orange band on illumination (光照) by UV light.


EB用量0.5 μg/mL , 50 ng的微量DNA可成带



> Application

- 1 Isolation of nucleic acids
- 2 Recycling DNA fragment
- 3 Estimating the approximate molecular weight of nucleic acids
- 4 Blotting





Gel excision and purification

SanPrep柱式DNA胶回收试剂盒


磁珠法DNA胶回收试剂盒



2.2.2 DNA libraries

- DNA libraries are sets of DNA clones, each of which has been derived from the insert of a different fragment into a vector followed by propagation in the host.
 DNA文库是一套DNA克隆,每个克隆都是插入了不同片段的载体在宿主中扩增后的产物。
- A clone Propagation of the host cell containing the recombinant DNA forms a set of genetically identical organism.

单克隆:带有某一种重组DNA的宿主细胞,增殖构成的一群具有遗传一致性的个体。



(1) Genomic libraries (基因组文库)

- Genomic libraries are prepared from random fragments of genomic DNA.
- Physical shearing and restriction digestion;
- The approach may be an inefficient method of finding a gene, especially in eukaryotic genomes, where much of the DNA is noncoding.



(2) cDNA library

- \rightarrow mRNA \rightarrow cDNA \rightarrow dsDNA \rightarrow Vector \rightarrow Host
- cDNA library are derived from the mRNA by reverse transcription and are then inserted into a vector.
- cDNA libraries are efficient for finding and cloning a gene, but only the coding region, not the surrounding genomic sequences.

筛选文库(screening libraries):寻找目的基因




Nucleic acid hybridization - Southern

blotting

A radioactively or fluorescently labeled DNA probe

Identification of the protein product

Antibody or enzyme activity measurement





2.3 DNA recombination

2.3.1 Restriction digestion (酶切)

将载体和扩增后的目的基因片段用相应的 限制性内切酶消化—电泳,回收酶切产物

2.3.2 Ligation


- (1) 具粘性末端的连接√
- (2) 平末端的连接
- (3) <u>TA克隆</u>
- (4) 同聚物加尾后的连接
- (5) 衔接物或人工接头的连接


-(DNA文库)





单酶切的连接问题

- Religation of the original vector plasmid.
 载体自身环化
- The target DNA fragment can be inserted into the vector with either orientation (方向) to form a recombinant molecule.











> 双酶切定向连接

Double digestion


同时解决载体自连








TA克隆


普通的Taq酶可在产物的3′端多加一个A,这是非模板依赖的。












衔接物的连接

符接物:是指用化学方法合成的一段由10-12 bp组成的、具有一个或数个限制酶识别位点的寡核苷酸片段。




.5'CCGAATTCGG

3'GGCTTAAGCCS'

EcoR I 核酸内切限制酶

5'AATTCGGCG3'

3'GCC CHARLES CONTRACTOR CONTRACT



人工接头的连接

OHGATCCCCGGG GGGCCC

DNA ligase


直接具备某一限制性内切酶 黏性末端的寡 核苷酸片段。

5'-OHG AT C C C C G G G G G G C C C C T A G_{OH}-5'



亚克隆 (Subcloning): 从载体到载体

Subcloning is the process to transfer of a fragment of cloned DNA from one vector to another.



- Express the gene in a particular species
- Enables us to investigate a short region of a large cloned fragment in more detail.



2.4 Introduce recombinant DNA into host cells

2.4.1 Transformation (转化)

Transformation is a process of uptake of exogenous DNA (normally plasmids) by a host (bacteria or yeast) cell.

(1) Competent cells (感受态细胞)

Competent cells are specially treated cells that foreign DNA can enter easily.

Calcium ions (Ca²⁺) treatment

Electrical treatment




(2) Transformation process

大肠杆菌(E. coli)转化过程:


重组DNA与competent cells在 冰上混合 (Taking up DNA)→ 42 °C heat-shock (热激) 1~2 min →加适量液体培养基培养, 平板(含抗生素)选择培养

Heat-shock induces enzymes involved in the repair of DNA, which allow the cells to recover from unusual conditions of the transformation process, and increases the efficiency.





Selection of plates (平板选择培养)





(3) Transformation efficiency (转化率)

Transformation efficiency: The number of colonies (菌落数) formed on a selective plate per µg transformed DNA (pure plasmid).





- > 10⁵ colonies/μg plasmids: would be enough for a simple cloning experiment
- >>10° colonies/µg plasmids: can be used for generation of libraries

(4) Factors that influence transformation efficiency

- Concentration, purity and conformation of transformed DNA 浓度越大、纯度越高,转化率越高;超螺旋>环状 >线型,分子大的转化率比分子小的要低。
- ➤ Physiological state of transformed cells CaCl₂处理贮藏之后的成活率。对数生长期。
- Environmental conditions of transformation 如温度,pH,离子浓度等。
- Restriction modification system of receptor cells
 - 一般用作转化受体的大肠杆菌菌株都是选用在宿主细胞控制的限制体系上有缺陷的突变体。



2.4.2 Transfection (转染)

- The introduction of DNA into mammalian cells is called transfection.
 外源DNA导入到哺乳动物细胞的过程称为转染(相当于细菌的转化)。
- Bacteria can be transformed with DNA extracted from a bacterial virus, a process known as transfection.
 外源DNA通过病毒载体的感染转移到宿主细胞的过程称为转染。

Transduction (转导): 通过病毒将一个宿主的DNA转移到另一个宿主的细胞中而引起的基因重组现象。



2.5 Screening and analysis of recombinants

Screening transformants: Allows all cells to grow, but tests the resulting clones for the presence of the insert DNA in the vector.

Nucleic acid hybridization (screening libraries)

Selectable markers

PCR


Ristriction digestion and electrophoresis

Sequencing




2.5.1 Blue-white screening (蓝白斑筛选)


- Insertional inactivation of the lacZ'
- Plate containing IPTG and X-gal
- The X-gal is converted to a blue product if the lacZ' gene is intact (完 整的) and induced by IPTG.
- Recombinants grow as white colonies.



2.5.2 Ristriction digestion and electrophoresis







- 1. Plasmid digest by EcoRI and Notl
- 2. Plasmid
- 3~5. Plasmid with exogenous DNA digest by *Eco*Rl and *Not*l




2.5.3 Sequencing

—TACGddA

-TddA






—TACGAddC

—TAddC

利用 ddNTP

作为DNA合 有特定末端的 DNA片段。




-TACddG

dNTP: ddNTP 的比例最好为

−ddT

1:100



Proliferation and preservation of transformants

- Grow in liquid broth (液体培养基) with antibiotics.
- Store a portion of culture with 15-30% glycerol (甘油) at -20℃ or -70℃.

The effect of glycerol stock is to protect the cells from ice crystal formation.




2.6 Expression of cloned genes

- cDNA克隆
- 可诱导的启动子(IPTG诱导)
- SD/Kozak序列
- 报告基因 (reporter gene): 编码一种在目标细胞中不正常产生的但容易被检测和鉴定的蛋白质基因。

检测方法:

- 酶活检测(酶蛋白)
- ➣ 免疫化学检测法 (Western blot)






SDS-PAGE





转膜



Western Blot 通 过SDS-PAGE区分 标记的二抗 显色(显影)对 靶蛋白含量进行



Summary

- 1. Definition and general protocol of DNA cloning
- 2. Types and properties of vectors
- 3. Types and properties of restriction endonucleases
- 4. Directional ligation of target DNA and vector
- 5. Transformation and transfection
- 6. Methods for screening recombinants







怎样将一个平末端DNA片段插入到载体 EcoR I 限制位点中去?